

Gel Extraction of DNA



E.Z.N.A.™ Gel Extraction Kit

1. Excise gel slice containing DNA fragment of interest.

- a. Gel electrophoresis fractionates DNA fragments.
- b. The gel is exposed to UV to find the DNA fragments (stained by Ethidium bromide).
- c. The goal DNA band is identified.
- d. Physically remove the slice of gel contains the goal DNA with clean surgical blade.

2. DNA Purification

- e. Determine the appropriate volume of the gel slice by weighing it in a clean 1.5 ml microcentrifuge tube.
- f. Add an equal volume of Binding Buffer (XP2).
- g. Incubate the mixture at 55°C-60°C for 7 min or until the gel has completely melted.
- h. Mix by shaking or vortexing the tube in increments of 2-3 minutes.
(IMPORTANT: Monitor the pH of the Gel/Binding Buffer mixture after the gel has completely dissolved. DNA yields will significantly decrease when pH >8.0. If the color of the mixture becomes orange or red, add 5 µl of 5M Sodium Acetate, pH 5.2 to bring the pH down. After this adjustment, the color of the Gel/Binding Buffer mixture should be light yellow.)
- i. Place a HiBind® DNA column in a provided 2 ml collection tube.
- j. Apply 700 µl of the DNA/agarose solution to the HiBind® DNA column, and centrifuge at 10,000 x g for 1 min at room temperature.
- k. Discard liquid and place the HiBind® DNA column back into the same collection tube. For volumes greater than 700 µl, load the column and centrifuge successively, 700 µl at a time. Each HiBind® DNA column has a total capacity of 25 µg DNA. If the expected yield is larger, divide the sample into the appropriate number of columns.
- l. Add 300 µl of Binding Buffer (XP2) into the HiBind® DNA column. Centrifuge at 10,000 x g for 1 min at room temperature to wash the column. Discard the flow-through and re-use the collection tube.
- m. Wash the HiBind® DNA column by adding 700µl of SPW Wash Buffer diluted with absolute ethanol. Centrifuge at 10,000 x g for 1 min at room temp.
(Note: SPW Wash Buffer Concentrate must be diluted with absolute ethanol before use. See label for directions. If refrigerated, SPW Wash Buffer must be brought to room temperature before use.)
- n. Repeat step 8 with another 700µl of SPW Wash Buffer diluted with absolute ethanol.

- o. Discard liquid and centrifuge the empty HiBind® DNA column for 2 min at maximal speed (13,000 x g) to dry the column matrix. Do not skip this step, it is critical for the removal of ethanol from the HiBind® DNA column.
- p. Place a HiBind® DNA column into a clean 1.5 ml microcentrifuge tube. Add 30-50µl (depending on desired concentration of final product) of Elution Buffer (10 mM Tris-HCl, pH 8.5) directly onto the column matrix and incubate at room temperature for 1 minute. Centrifuge for 1 min at maximal speed (13,000 x g) to elute DNA. This represents approximately 70% of bound DNA. An optional second elution will yield any residual DNA, though at a lower concentration.