

In-Fusion Cloning

RECOMMENDED IN-FUSION REACTIONS FOR PURIFIED FRAGMENTS			
Reaction Component	Cloning Reaction	Negative Control Reaction	Positive Control Reaction
Purified PCR Fragments	10-200 ng*	-	2 ul of 2 kb control insert
Linearized Vector	50-200 ng**	1 ul	1 ul of pUC19 control vector
5x In-Fusion HD Enzyme Premix	2 ul	2 ul	2 ul
Deionized water	To 10 ul	To 10 ul	To 10 ul

* <0.5 kb: 10-50 ng; 0.5 to 10 kb: 50-100 ng; >10 kb: 50-200 ng
** <10 kb: 50-100 ng; >10 kb: 50-200 ng

1. Set up the In-Fusion cloning reaction:

5x In-Fusion HD Enzyme Premix	2 ul
Linearized Vector	x ul*
Purified PCR Fragment	x ul*
dH ₂ O (as needed)	x ul

Total Volume 10 ul

** For reactions with larger volumes of vector and PCR insert (> 7ul of vector + insert), double the amount of enzyme premix, and add dH₂O for a total volume of 20ul.*

- Adjust the total reaction volume to 10 ul using deionized H₂O and mix the reaction.
- Incubate the reaction for **15 minutes at 50°C**, then place on ice.
- Continue to the transformation procedure. You can store the cloning reactions at -20°C until you are ready. You can use 2-5ul of the reaction to transform.