

## Restriction Digests

Use 50 $\mu$ l reactions for plasmid digestion. Use the smaller reaction for diagnostic purposes.

- 5 $\mu$ l buffer #X (check with NEB)

- 0.5 $\mu$ l BSA (if you need it)

- 35-44 $\mu$ l PCR product

- water up to 50 $\mu$ l

- 1 $\mu$ l enzyme

= 50 $\mu$ l

- 2.5 $\mu$ l buffer #X (check with NEB)

- 0.25 $\mu$ l BSA (if you need it)

- 17-22 $\mu$ l PCR product

- water up to 25 $\mu$ l

- 0.5 $\mu$ l enzyme

= 25 $\mu$ l

1. Check for specific temperatures and times for optimal digestion of your enzyme (neb.com).
2. Vortex, spin down
3. 37 °C incubator, 30 to 60 min.
4. Heat inactivate.

### TIPS:

- Always keep enzyme in -20 °C bucket.
- When doing multiple RDs: make a Master Mix (everything except DNA), vortex, divide into your reactions, and finally add DNA.