

Gibson Assembly

For assembly, the following was found to be highly reliable

1. Prepare the Assembly Mix on ice
2. Incubate for one hour at 50°C
3. Transform

Assembly Mix	Volume (μ l)
Gibson Master Mix	15
Total DNA	5
Total	20

Master Mix recipe

Gibson Master Mix	Volume (μ l)
Taq ligase (40u/ μ l)	50
5x isothermal buffer	100
T5 exonuclease (1u/ μ l)	2
Phusion polymerase (2u/ μ l)	6.25
Nuclease-free water	216.75
Total	375

Whilst the mixture can be refrozen, for best performance you should aliquot it into 25x15 μ l.

5x isothermal buffer	Volume (μ l)
25% PEG-8000	0.75g
500 mM Tris-HCl pH 7.5	1500
50mM MgCl ₂	75
50mM DTT	150
1mM dATP	30
1mM dTTP	30
1mM dCTP	30
1mM dGTP	30
5mM NAD	300
Nuclease-free water	...
Total	3000

Use nuclease-free water to make up to 3000 μ l. In this case, very little water is needed as the volume is 2145 μ l + 0.75g.