

## **Protocols for Transformation of bacteria via TSS Method**

iGEM Grenoble-EMSE-LSU 2013  
Stewart Humble ([stewhum1@me.com](mailto:stewhum1@me.com))  
June 2013

### ***Adapted for cultures in Falcon (15 mL) Tubes***

- 1. Make a 1mL bacteria pre-culture the day before, until saturation, in a 15mL Falcon Tube**
- 2. Next morning, dilute this culture 1/100 in a new 15mL Falcon Tube**
- 3. Incubate at 37°C in a shaker with the lids of the tubes loose for oxygen penetration into the culture**
- 4. Wait around 2-3 hours until OD600 reaches  $\approx 0.4-0.5$**
- 5. Prepare 0.5mL eppendorfs (in case heat shock is done with PCR machine) or 1.5mL eppendorfs (if heat shock is done using a 42°C water bath)**
  - *Put the eppendorfs on ICE*
  - *For each transformation:*
    - Pipette 50 $\mu$ L of cell culture per eppendorf
    - Pipette 50 $\mu$ L of TSS 2x per eppendorf; Mix WELL but GENTLY
- 6. Incubate on ice for 45 minutes**
- 7. Add DNA into each eppendorf according to each transformation (in case of a biobrick, add 1-2 $\mu$ L of the red DNA solution)**
- 8. Incubate on ice for 10 minutes**
- 9. Heat shock at 42°C by using the PCR machine or water bath for 90 seconds**
- 10. Incubate on ice for 10 minutes**
- 11. Add 900 $\mu$ L of LB (0.2% glucose is optional) for each eppendorf (if using 0.5mL tubes then transfer solution to new COLD 1.5mL tubes)**
- 12. Incubate at 37°C for 1-2 hours**
- 13. Plate on LB-Agar supplemented with the appropriate antibiotic. Wait until the next day to harvest colonies.**