

Serum Incubation of Minicells

Materials:

- **Human Serum from human male AB plasma, USA origin, sterile filtered**
- **Phosphate-buffered saline (PBS) 1X solution**
1 L volume, store at room temperature
 - NaCl (8 g)
 - KCl (0.2 g)
 - Na₂HPO₄ (1.44 g)
 - KH₂PO₄ (0.24 g)Dissolve reagents in 800 mL of H₂O, adjust pH to 7.4 with HCl, and then add H₂O up to 1L.
- **SDS sample buffer**
30 mL volume, store at 4°C
 - 1M Tris pH 6.8 (0.72 mL)
 - 100% glycerol (3 mL)
 - 10% SDS (6 mL)
 - β-mercaptoethanol (1.5 mL)
 - 1% bromophenol blue (3 mL)
 - ddH₂O (15.78 mL)

Directions:

1. When the plasma ships in, freeze in small aliquots (5 mL each) in the -80 freezer. Don't freeze-thaw these too many times, or they will lose activity.
2. Setting up the control: Heat inactivate a small portion (one or two samples, where each sample ~500 uL) of serum for ~30 minutes at 56°C.
3. Incubate minicell with % serum 0-100% (by volume with buffer = phosphate buffered saline (PBS)) for 20 minutes at 37°C, mixing every 5 minutes.
4. Pellet minicells with Kay's centrifuge at 10,000xg. Keep supernatant as control- run gel immediately- freeze if gel won't be run immediately.
5. Wash minicells in PBS twice- 5 minutes each- let it resuspend, pellet again, resuspend pellet in PBS, add equal volume of 2x SDS sample buffer and boil for 5 minutes, freeze, then boil again.
6. Before loading gel, thaw samples, boil again for 3 minute, microfuge for 1 minute, take sample out from top Run an SDS page gel with solubilized washed minicells. Cut gel in half: stain one half for visualization with Coomassie blue, perform a western blot on the other half with C3 antibodies to check for complement activation.

Western Blot (standard)

Transferring to Nitrocellulose Membrane:

Materials:

- Transfer Buffer- Tris-glycine buffer (stored in -20)
 - 10% 10X Tris-glycine stock solution
 - 20% Methanol

Serum Tests for Complement Deposition

- 70% water

Need about 1.4 L per box (that means 140 mL Tris-glycine, 280 mL of methanol, and 980 mL water)

Directions:

1. Soak membrane in methanol for 15 seconds
2. Rehydrate membrane in water for 1 minute
3. Soak membrane in transfer buffer for 5 minutes
4. Soak Blotting Paper(s) in transfer buffer
5. Cut out gel
6. Assemble sandwich in the following order
 - a. Clear plate
 - b. Sponge
 - c. Blot paper (roll out bubbles after placement)
 - d. Wet gel with transfer buffer *prevents tearing* (roll out bubbles after placement)
 - e. Nitrocellulose membrane (roll out bubbles after placement)
 - f. Second piece of blot paper (roll out bubbles after placement)
 - g. Sponge
 - h. Clear Plate
7. Close sandwich and add ice
8. Run at 200 milliamps for 2 hrs. or 30-100 milliamps overnight
Runs to red(+)

Blocking:

Materials:

- Blocking Buffer- 5% non fat milk in 1xPBS at a 1:1 ratio (you need at least 50 mL to cover the whole thing)
- PBST- 1X PBS + 0.05% Tween: recommended to make about 500 mL (500 mL PBS with 250uL Tween)

Directions:

1. Wash 2-3 times in deionized water and mark lanes and markers with a ballpoint pen
2. Typically block overnight at 4 degree or 2 hrs. at room temperature
3. Dilute primary antibody solution with blocking buffer to appropriate concentration
**use 1:1000 dilution*
4. Add 50 mL to membrane- just enough to cover it (or do it in a bag to save antibodies)
5. Leave in the 4 degree for 2-3 hours
 - a. If in fridge, take out every 30 minutes and shake back and forth
 - b. If in cold room, leave on bellydancer
6. Wash twice quickly, twice for 5 minutes and twice for 10 minutes
7. Dilute secondary antibody with blocking buffer
**use 1:5000 dilution*
8. Incubate in room temperature (on belly dancer if possible) for 30 min- 1 hr (longer incubation time leads to longer band visibility)
9. Wash 4 times, 5-10 minutes each with PBST
10. Wash 2 times with PBS for 5-10 minutes