

6-18-13

Plasmid Prep of J04450, K299009, K176012, K412005
and J04450(unknown 6)

* Uses 70% ethanol to precipitate the DNA instead of
100% ethanol.

- Put in Freezer for 30 minutes, and spin in ~~microcent~~
microcenter fuge for 5 minutes.

NB, CL, GE

6-18-13

Finishing the Cadmium tests from page 12.

- Spin down tubes
- Pour off liquid into beaker.
- Add 1 mL of 1x PBS to each tube and resuspend pellet.
- Transfer liquid to cuvettes.

Sample	M[Ca]	OD ₆₀₀	Emission 509 Fluorescence 405	Emission 380.1 Fluorescence 405 (Dilution)
1	0	1.088	0.049	0.689 (1 to 10)
2	.001	1.199	0.054	0.781 (1 to 10)
3	.005	0.397	0.031	0.243 (1 to 10)
4	.01	0.585	0.037	0.400 (1 to 10)
5	.05	1.095	0.069	0.092 (1 to 100)
6	.1	1.051	0.038	0.150 (1 to 100)

	Fluor/OD
1	6.333
2	6.514
3	6.121
4	6.838
5	8.402
6	14.272

6-18-13

- Started Cadmium Experiment with liquid culture from 6/11.
- Added 4 ml of LB + chlor to both positives and negatives (4 λ of chlor)
 - Added 100 λ of Cadmium Detector cells.
 - Added Cadmium according to chart.

Tubes	1	2	3	4	5	6
Molarity	0	.001	.005	.01	.05	.1
Cadmium (ul)	0	4	20	40	200	400

- Incubate

* Started a new culture

	Sample m(c)	OD ₆₀₀	Emission ³⁶⁰⁰¹ Fluorescence 405	Fluor/OD
1	0	774.700	0.980 (1 to 10)	2.5
2	.001	.780	0.226 (1 to 10)	3.013333333333
3	.005	.067	0.175	2.61194029851
4	.01	.256	0.187 (1 to 10)	7.3046875
5	.05	1.420	0.074 (1 to 100)	5.21126760563
6	.1	1.611	0.151 (1 to 100)	9.37306021105