


## USER Cloning

By Penn iGEM 2013 and Spencer Glantz


**Goal:** *To scarlessly combine two linear fragments of DNA, each of which contain uracil residues, into a larger linear or circular fragment.*

See figures below that explain the concept behind USER cloning (for a hypothetical situation where one wants to scarlessly ligate Part1 to Part2)


### PCR **Part 1** off of Template 1



## PCR Part 2 off of Template 2



## Digest PCR Products with DpnI and User Enzyme



## Protocol:

(1) Amplify target sequences (Parts 1 and 2) with primers that contain an adenine at the 5' end and a uracil roughly 8-11 bp later (in place of a thymine). In between the 5' adenine and the uracil will be the USER cloning overhang sequence - which may be standardized or part of the gene you are trying to clone. When designing overhangs, I try to find an A and T roughly 8-11 bp away from each other with a lower GC content (higher specificity). Use Pfu Turbo Cx polymerase to minimize PCR errors - although I have only had success with this polymerase for amplicons  $\leq$  2 kb.

### Pfu Turbo Cx Amplification Reaction Mix (20 $\mu$ L)

1  $\mu$ L Pfu Turbo Cx  
1  $\mu$ L 10  $\mu$ M Fwd Primer  
1  $\mu$ L 10  $\mu$ M Rev Primer  
1  $\mu$ L DNA template  
1  $\mu$ L 10 mM dNTPs  
2  $\mu$ L Pfu Turbo Cx Buffer (10x)  
13  $\mu$ L H<sub>2</sub>O

### Thermocycler Settings:

95°C for 3 minutes (or 5 minutes for colony/genomic DNA templates)  
30 Cycles of:  
    95°C for 30 seconds  
    55°C for 30 seconds  
    68°C for extension time (1 kb/minute) + some (I'm generous here)  
1 Cycle of:  
    68°C for 20 minutes  
    Hold at 4°C

(2) Perform USER cloning reaction:

\*Note, run reactions with just each part on its own as negative controls to see how much "background" you are getting from each one - replace missing parts with H<sub>2</sub>O. This is most important for parts that contain an origin of replication and an antibiotic resistance marker\*

2-way assembly Reaction Mixture (4.8 uL)

0.4 uL DpnI  
0.4 uL USER enzyme  
2 uL Part 1  
2 uL Part 2

3-way assembly Reaction Mixture (4.8 uL)

0.4 uL DpnI  
0.4 uL USER enzyme  
1.33 uL Part 1  
1.33 uL Part 2  
1.33 uL Part 3

On the thermocycler:

37°C for 1 hour.  
22°C (room temperature) for 15 minutes.  
Cool to 4°C.

Directly transform 1 uL of the product into chemically competent cells (follow plasmid transformation protocol) - Plate everything and should get 20 - 50 colonies.