

Colony PCR

Preparing for *E. coli* colony PCR:

1. Resuspend colony or pellet from overnight culture in 10-15ul ddH₂O.
2. Set up a PCR reaction:

Water	to a total of 25ul
5X Buffer*	5.0ul
10mM dNTPs	0.5ul
Template	5.0ul
5' Primer	0.5ul
3' Primer	0.5ul
Phusion Polymerase	0.1 ul
	= 25ul reaction

* DMSO can be added for troublesome PCR's (~0.1ul).

3. Use 5ul of the resuspended colony as the template for the PCR.
4. Place tubes in thermocycler and set the times according to this table (3-step protocol):

Cycle step	2-step protocol		3-step protocol		Cycles
	Temp.	Time	Temp.	Time	
Initial denaturation	98°C	30 s	98°C	30 s	1
Denaturation	98°C	5-10 s	98°C	5-10 s	25-35
Annealing*	-	-	X°C	10-30 s	
Extension	72°C	15-30 s/kb	72°C	15-30 s/kb	
Final extension	72°C	5-10 min	72°C	5-10 min	1
	4°C	hold	4°C	hold	

** IMPORTANT: For the **initial denaturation** step, the time should be 3-5 minutes in order to break open the cells. Everything else should be the same. **

5. After PCR has been completed, add dye to samples and load them on a gel to see if the insert is there.