


Colony PCR

1. Keep everything on ice. The PCR reaction mixture in the ice preparation, and then placed in a PCR reaction in the PCR reaction apparatus. This cold start (Cool Start Method) can enhance the specificity of PCR amplification to reduce non-specific reactions in the PCR process, can get good PCR results.
2. Make up a master mix of everything into one microcentrifuge tube.
3. Pipette up and down in the microcentrifuge tube, drain 25 μ L or 50.0 μ L solution to each PCR tube.

a. **25.0 μ L reaction system**

18.3 μ L ddH₂O
2.5 μ L 10x Buffer (mg²⁺)
2.0 μ L dNTPs
0.5 μ L forward primer
0.5 μ L reverse primer
0.2 μ L DNA polymerase
Colony stab (template DNA)

-----**25.0 μ L** Total



b. **50.0 μ L reaction system**

37.75 μ L ddH₂O
5.0 μ L 10x Buffer (mg²⁺)
4.0 μ L dNTPs
1.0 μ L forward primer
1.0 μ L reverse primer
0.25 μ L DNA polymerase
Colony stab (template DNA)

-----**50.0 μ L** Total

4. Pick colonies from plates, spot onto these PCR tubes.
5. Run the "Colony PCR" program, and adjust your extension time as described.

The "Colony PCR" program

Initial denaturation: 95°C for 5:00min

25 cycles of:

94°C for 0:30 min

55°C for 0:30 min (different primers different annealing temperature)

72°C for t min ("t" depends on the length of goal sequence, 1min per 1kb)

Final extension: 72°C for 10:00 min