


# Colony PCR

1. Keep everything on ice. The PCR reaction mixture in the ice preparation, and then placed in a PCR reaction in the PCR reaction apparatus. This cold start (Cool Start Method) can enhance the specificity of PCR amplification to reduce non-specific reactions in the PCR process, can get good PCR results.
2. Make up a master mix of everything into one microcentrifuge tube.
3. Pipette up and down in the microcentrifuge tube, drain 25 $\mu$ L or 50.0 $\mu$ L solution to each PCR tube.

a. **25.0 $\mu$ L reaction system**

**18.3 $\mu$ L** ddH<sub>2</sub>O  
**2.5 $\mu$ L** 10x Buffer (mg<sup>2+</sup>)  
**2.0 $\mu$ L** dNTPs  
**0.5 $\mu$ L** forward primer  
**0.5 $\mu$ L** reverse primer  
**0.2 $\mu$ L** DNA polymerase  
Colony stab (template DNA)

-----**25.0 $\mu$ L** Total



b. **50.0 $\mu$ L reaction system**

**37.75 $\mu$ L** ddH<sub>2</sub>O  
**5.0 $\mu$ L** 10x Buffer (mg<sup>2+</sup>)  
**4.0 $\mu$ L** dNTPs  
**1.0 $\mu$ L** forward primer  
**1.0 $\mu$ L** reverse primer  
**0.25 $\mu$ L** DNA polymerase  
Colony stab (template DNA)

-----**50.0 $\mu$ L** Total

4. Pick colonies from plates, spot onto these PCR tubes.
5. Run the "Colony PCR" program, and adjust your extension time as described.

**The "Colony PCR" program**

Initial denaturation: 95°C for 5:00min

25 cycles of:

94°C for 0:30 min

55°C for 0:30 min (different primers different annealing temperature)

72°C for t min ("t" depends on the length of goal sequence, 1min per 1kb)

Final extension: 72°C for 10:00 min