

## Ligation Reactions (20 $\mu$ L)

1. Gel purify all digested inserts and plasmid backbones
  - o Refer to gel purification of DNA protocol
2. Add the following reaction components to a 0.2 mL PCR tube
  - o 11  $\mu$ L - ddH<sub>2</sub>O
  - o 2  $\mu$ L - upstream insert
  - o 2  $\mu$ L - downstream insert
  - o 2  $\mu$ L - destination plasmid (backbone)
  - o 2  $\mu$ L - T4 DNA Ligase Buffer (10x)
  - o 1  $\mu$ L - T4 DNA Ligase
3. Incubate in thermocycler w/ heated lid
  - o Overnight @ 16° C
  - o Heat inactivate for 20 min. @ 80° C (optional)

### Alternative protocol

1. Gel purify all digested inserts and plasmid backbones
  - o Refer to gel purification of DNA protocol
2. Add the following reaction components to a 0.2 mL PCR tube
  - o 5  $\mu$ L - ddH<sub>2</sub>O
  - o 4  $\mu$ L - upstream insert
  - o 4  $\mu$ L - downstream insert
  - o 4  $\mu$ L - destination plasmid (backbone)
  - o 2  $\mu$ L - T4 DNA Ligase Buffer (10x)
  - o 1  $\mu$ L - T4 DNA Ligase
3. Incubate for ~10 min. @ room temperature
4. Heat inactivate for 20 min. @ 80° C (optional)